204

Manual Examination Table

Ritter®

[English](#)

[Deutsch](#)

[Français](#)

[Español](#)

[Dansk](#)

[Suomi](#)

[Italiano](#)

[Svenska](#)

[Norsk](#)

[Nederlands](#)

[Polski](#)

[Türkçe](#)

ArtFull

User Guide

CE

0 0 3 - 1 0 3 6 3 - 9 9

Style U

003-10363-99 Rev AA7 (12/09/25)

TP202 20-42-FO-00014 Rev A1 C2169

Product Information

(The information below is required when calling for service)

Dealer [name / phone]:

Date of Purchase:

Model / Serial Number:

***Midmark Authorized Service
Company*** [name / phone]:

***Model / Serial
Number Location***

ArtFull

Table of Contents

Important Information

- Safety Symbols..... 4
- Warnings..... 5
- Intended Use 5
- Electrical Requirements..... 5
- Electromagnetic Interference..... 6
- Disposal of Equipment..... 6
- Transportation / Storage Conditions 6
- Authorized Representatives 6

Operation

- Back Section Positioning 8
- Stirrups 9
- Footrest & Treatment Pan..... 10
- Footstep..... 10
- Paper Roll & Tear Strap..... 11
- Pelvic Tilt (optional) 12
- Receptacle (optional)..... 13

Maintenance

- Cleaning 14
- Periodic Maintenance 14
- Calling For Service 14

Specifications

- Weights & Dimensions..... 15
- Upholstery Dimensions..... 16
- Voltage, Current, Frequency Ratings & Compliance Chart 17

Warranty Information

- Scope of Warranty..... 18

Important Information

Safety Symbols

WARNING
Indicates a potentially hazardous situation which could result in serious injury.

Caution
Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert
Indicates a potentially hazardous situation which could result in equipment damage.

Note

Amplifies a procedure, practice or condition.

Proper shipping orientation

Fragile

Type B, Applied Part
(Upholstery)

Protective earth ground

Patient Weight Limit

Manufacturer


Serial Number

Keep dry

Maximum stacking height
(**palletted** units)

Fuse rating specification

Dangerous voltage / shock hazard

Refer to manual

Catalogue Number

Medical Device

Warnings

WARNING

To avoid injury, do not put head or extremities under the paper tear strap.

WARNING

Connecting equipment to the multiple socket outlet effectively leads to creating a Medical Electrical System and the result can be a reduced level of safety.

WARNING

Table must be positioned in normal use such that the power cords can be easily accessed to unplug.

WARNING

To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

WARNING

This product can expose you to chemicals including Antimony Oxide (Antimony Trioxide), which is known to the state of California to cause cancer. For more information go to www.p65warnings.ca.gov

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

WARNING

No modification of this equipment is allowed.

Intended Use

This product is intended to be used as a table to provide positioning and support of patients during examination procedures conducted by qualified medical professionals in an office environment.

Electrical Requirements

Caution

To completely isolate the table from electrical mains supply, both power cords must be unplugged.

When using high frequency surgical devices or endocardial catheters:

- Use non-conductive material to insulate patient from metal portions of the table.*
- Consult operating instructions for the device before using in conjunction with drawer or upholstery heater.*

Failure to comply may result in electric shock or burns to patient.

Electromagnetic Interference

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room.
- Plug table into an isolated circuit.
- Increase separation between table and interfering device.
- Contact Midmark if interference persists.

Disposal of Equipment

At the end of product life, the table, accessories, and other consumable goods may become contaminated from normal use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.

Transportation / Storage Conditions

Ambient Temperature Range:.....-5°C to +38°C (+23°F to 100°F)
Relative Humidity.....10% to 90% (non-condensing)
Atmospheric Pressure500hPa to 1060hPa (0.49atm to 1.05atm)

Authorized Representatives

EU

REP

WARNING

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in the EU should direct all questions, incidents, and complaints to Midmark’s Authorized Representative listed below.

AR Experts BV
Boeingavenue 209
1119 PD Schiphol-Rijk
The Netherlands
www.ar-experts.eu

Customers in the UK should direct all questions, incidents, and complaints to Midmark’s UK Responsible Person listed below.

UKCA Experts Ltd
6th Fl. City Gate East,
Toll House Hill,
Nottingham, NG1 5FS,
UK

Authorized Representatives - continued

Customers in India should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Midmark India

Kalpataru Point, Unit No. 12, 1st Floor Kamani road, Sion (E)

Mumbai, 400022

Tel: +91 22 4915 3000

Fax: +91 22 4915 3100

Customers in Saudi Arabia should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Attieh Medico Ltd

Abdullah M. Al Khalifi Street

Al Nakheel Dist. II

P.O. Box 116105

Jeddah (East)

21391

Saudi Arabia

Tel: + 966 2 286 4707

Fax: + 966 2 286 4744

Customers in the United Arab Emirates should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

GULF DRUG LLC.

Al Barsha,

Dubai, UAE

Phone: +971 4 501 4000

Fax: +971 4 501 4100

Customers in Qatar should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Tadmur Trading

Midmac Roundabout

Opposite Al Mana Petrol Station

Doha, Qatar

Tel: +974-44337000

Fax: +974-44337100

Customers in Israel should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Medtechnica Ltd.

HaTnufa St 7

Petah Tikva, Israel

Tel: +972 3-925-4040


Fax: +972 3-924-9977

Operation


Back Section Positioning

To position the back section...

- A) While holding the back section in place, squeeze the Back Release Handle (either side).
- B) Move the back section to desired position.
- C) Release handle to lock back section in position.

Back Release Handle

MA697500i

Stirrups


Caution

*Be sure the stirrups are locked in place before using.
The stirrups will not support the patient's entire weight.*

To position the stirrups...

- A) Pull the stirrup out, then unfold.*
- B) Lift the stirrup slightly, then move it left or right as desired.*
- C) Release stirrup to lock in position.*

Footrest & Treatment Pan

Caution

Do not use footrest to support patient's entire weight.

To access the treatment pan...

Lift footrest pad slightly, then slide pad back into table.

To extend footrest...

Pull footrest out to desired position.

Footstep


Caution

Ensure patients place their feet fully on the flat surface of the footstep when getting aboard and exiting the table.

To access the footstep...

Lift footstep slightly, then pull outward.

Paper Roll & Tear Strap

Note	
Paper Roll Size (max.):	
Seamless Upholstery.....	21 in. wide x 3,5 in. diameter (53 cm x 9 cm)
Soft Touch Upholstery	21 in. wide x 4.75 in. diameter (53 cm x 12 cm)

Pelvic Tilt (optional)

WARNING



To avoid risk of injury do not put extremities under the seat.

To elevate seat section...

Lift seat until pelvic support drops into place.

To return seat section to flat position...

- A) Lift seat section slightly.
- B) Pivot trigger handle upward.
- C) Lower seat section.



Receptacle (optional)

The duplex receptacle on the side of the table, provides power for accessories used during medical procedures. There are two circuit breakers located under the seat section. If the receptacle's maximum load is exceeded, the circuit breakers interrupt power to the receptacle.

Maximum Load 5 amps

To reset circuit breaker(s)...

Push the circuit breaker toggle switch.

MA9894p

Maintenance

Cleaning

Upholstery

Equipment Alert

The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery.

Wash your upholstery weekly with a mild liquid soap and water mixture, rinse with clear water and dry completely to remove disinfection cleaner build-up.

Disinfect your upholstery with a solution of standard bleach and water mixed 1 in 10 (10%) or chlorine based cleaners. Follow this with a clear water rinse and thorough drying of material. See current CDC Guideline for Disinfection & Sterilization in Healthcare Facilities.

To minimize damage caused by disinfectant cleaner residue build-up, do not allow disinfectants to pool on the upholstery surface. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

Detailed care and maintenance instructions are included with your product. This information is also available on [midmark.com](https://www.midmark.com) in the Technical library under the User Information tab for your product.

Painted Metal / Plastic Surfaces

Clean the painted metal and plastic surfaces weekly using a clean soft cloth, and mild cleaner.

Preventative Maintenance

Periodically inspect the following areas:

- Power cord should be free of cuts or other visible damage.
- All fasteners should be present and tightened securely.
- All mechanical functions should operate properly.

Have a qualified service technician inspect your equipment every twelve months.

Calling For Service

If service is required, contact your authorized Midmark dealer.

Note

Model / Serial Number information is required when calling for service.


If service is required, contact Midmark directly:

US 844.856.1230 / Canada 937.526.8585
Monday-Thursday 8:00 AM until 6:00 PM
Friday 8:00 AM until 5:00 PM (ET)
[midmark.com](https://www.midmark.com)

Specifications

Weights & Dimensions




Maximum Patient Weight:	500 lbs (226 kg)
Paper Roll (maximum size):	[See Paper Roll & Tear Strap page]
Range of Motion & Dimensions:	(See illustration below)
Weight of Table:	220 lbs (99.8 kg)
Power Cord:	6 ft. (1.8 m) long

MA8376i

Specifications (continued)

Upholstery Dimensions

MA83771

Specifications (continued)

Voltage, Current, Frequency Ratings & Compliance Chart

Model	Description				Electrical Ratings:		
		ES 60601-1	CAN / CSA C22.2, 60601-1	CE	VAC +/- 10%	Amps	Cycles (Hz)
204-005 204-013	Ritter® Exam Table - Export (Pass-thru Drawers) w/ Pelvic Tilt	•	•	•	-	-	-
204-006 204-014	Ritter® Exam Table - Export (Pass-thru Drawers) w/ Receptacle & Pelvic Tilt	•	•	•	230	5	50/60

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Midmark Corporation